
TECHNIQUES OF SYNTHESIS

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
² Phosphormolybdic acid (12 g) in ethanol (250 ml)

³ 4-Methoxybenzaldehyde (15 g) in ethanol (250 ml) + H₂SO₄ (2.5 ml)

⁴ 15 % aq. H₂SO₄ saturated with ceric sulphate

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Flash chromatography using dry packing and vacuum elution

VLC separation chromatography was performed on a column (D= 30/50 mm, h = 150 mm) with a suction joint in the bottom. The stationary phase was Merck Silica gel 60H. The conical bottom of the column was filled with cotton wool with a piece of column fitting filter paper on top. The column was loaded with SiO₂ under vacuum and the surface was pressed to uniform consistency and horizontal level by a cylindrical piece of glass. A piece of column fitting filter paper was placed on top of the SiO₂.

The sample (which should be weighed) was dissolved in the smallest possible amount of an appropriate solvent (typical DCM). ⁶Careful application on the column is crucial for obtaining a good resolution of fractions by the subsequent separation. The column was eluted with an appropriate solvent by applying suction to the column. (Optimum separation is obtained at R_f 0.3). However, gradient elution should be considered if fast running and slow running constituents have to be separated. The eluent can be applied using a reservoir on top of the column.

The fractions collected were checked by TLC. Fractions containing identical constituents were combined. Fractions were evaporated to dryness and weighed.

⁶ Alternatively, the crude product can be dissolved in an appropriate solvent, a small amount of SiO₂ added, the mixture evaporated to dryness and the residue applied in a horizontal layer on top of the SiO₂ in the column.

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
Screw-Cap Vessels for Medium Pressure Reactions and Processes Run under an Inert Atmosphere⁷

Reactions requiring temperatures higher than those of the boiling points of the solvents are traditionally run in expensive autoclaves or in sealed glass tubes. Sealing of glass tubes requires skill and caution. For small scale reactions at medium pressure (up to ca. 10 bar), inexpensive and read-to-make screw-cap closed reaction vessels are much easier to handle. Satisfactory results have been obtained with vessels made from commercial 13 mm glass screw-thread joints. Different types are presented in the figure and can be manufactured by any professional glassblower. The vessels are closed with standard Quick Fit 13-mm plastic screw caps equipped with metal and rubber disks. Butyl rubber and neoprene are most suitable, the latter being the more resistant to higher temperatures. In the case of aggressive solvents or reactants, the appropriate rubber type should be found by running a blank experiment. Alternatively, a disk of Teflon foil may be used to protect the rubber disk. The screw cap is carefully tightened with the finger tips. Using too much force in closing may lead to breakage of the glass screw thread.

For reasons of safety the vessels must always be protected or shielded during reactions at elevated temperatures or when hot. In the case of long, cylindrical vessels an iron water pipe closed at both ends with mufflers may serve. Reactions may be run in a thermostated heating cabinet or, when stirring is required, in a cylindrical aluminium block with borings for short cylindrical or pear-shaped vessels and a hole fitting a contact thermometer. The block, when furnished with an aluminium screw cap, may at the same time serve for protection. The contact thermometer is connected to the thermo-relay of a standard magnetic stirrer with heating, on which the aluminium block is placed. Proper current adjustment makes it possible to maintain reaction temperatures within the block to an accuracy of about ± 1 °C. Temperature equilibrium is reached slowly (after ca. 1 h) with pear shaped vessels placed in cylindrical borings unless a liquid heating transfer medium (e.g. paraffin oil) is placed in the boring. Reactions in methanol have routinely been run at 140 °C without any loss of material. Processes sensitive to moisture or oxygen can be run under an inert atmosphere in the screw-cap closed vessels using a metal disk with a hole. This allows perforation of the rubber disk with a syringe for evacuation as well as introduction of inert gas and reactants. With a 1-mm hole and a corresponding syringe, medium-pressure reactions can be run also after puncture of the system. In this case the use of two butyl rubber disks is recommended. If the septum has to be replaced under inert conditions, plastic bag techniques is used.

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Acid and Base properties

Acidity and Basicity in Non-activating solvents

Acids	pK _a	Bases	pK _b
ArSH	7		
CH ₃ NO ₂	9		
RSH	10	Et ₃ N	11
HOH	14	NaOH	16
		[1,1,1]Kryptand	16
ROH	15-17	NaOEt	18
CH ₃ COCl	16	1,8-(Me ₂ N) ₂ -Naphtalin	18
CH ₃ CHO	17	Pentamethylpiperidine (PMP)	19
		KOBu ^t	19
CH ₃ COMe	20		
C≡CH	25		
CH ₃ COOMe	25		
CH ₃ C≡N	25	Diazabicyclononen (DBN)	25
CH ₃ CONMe ₂		Diazabicycloundecen (DBU)	26
		LiN(TMS) ₂	26
		NaH	35
		NaNH ₂	35
		Li(Pr ⁱ) ₂	36
		LiTMP	37
PhH	41		

C=C-CH ₃	43		
C=CH ₂	44		
		Phospazen (tBu-P ₅)	45
C-CH ₃	50	LiBu ⁿ	50
		LiBu ^t	>50
		NaSiMe ₃	≈ 55

Usually a base providing a pK_a difference larger than 4 should be chosen by deprotonation.


Factors which influence deprotonation

Thermodynamic basicity	See table of basicities above
Kinetic basicity	Bu ⁿ Li and Bu ^t Li are slow MeLi, Bu ^s Li and PhLi are medium O- and N-bases are fast
Bulkyness	Consider Bu ^t Li, LiTMP, DBU, PMP, KOBu ^t
Reversible or non-reversible base	BuLi, RMgX, MetalH, CO ₃ ⁻² and SO ₃ ⁻² are irreversible
Charge	Uncharged bases are usually better soluble and kinetically faster
Polarity	Polarity decreases in the order Cs > K > Na > Li > Mg
Hardness or softness	Hardness decreases in the order O- > N- > C- > S
Nucleophilic reactivity	Consider Reaction with solvent Consider Reaction with electrophile (addition or substitution) Consider Reduction
Solubility	Consider solubility of base substrate and electrophile
Solvent	Consider different solvents for metallation and electrophilic addition in order to avoid side reactions
Accelerators	TMEDA for Bu ⁿ Li
Activation	LiH, NaH, KH: <ul style="list-style-type: none"> • Treatment with I₂/MeI • Washing with LiAlH₄ • Addition of crown-ether • Ultrasonification
Fresh preparation	LiH, NaH, KH
Base combinations	Base combinations can speed up deprotonation rate e.g. Bu ⁿ Li + KOBu ^t
Assistant base needed	Favorable properties may be combined e.g. strength, solubility, irreversibility

Factors which influence Reactivity and Regioselectivity by deprotonation followed by electrophilic addition

Nature of the base	Strong bases favor kinetic products Hard bases favor deprotonation for halogen-metal exchange Kinetic fast bases favor kinetic product Bulky bases favor thermodynamic product
Temperature	Low temperature favors kinetic product
Nature of electrophile	Hard electrophiles favor addition at hard sites
Order of addition	Reverse addition favors thermodynamic product
Nature of substituents	Consider directing properties Consider bulkyness

Strategy by deprotonation


Chemoselectivity and Regioselectivity of Organometallics M-R

Metal in M-R		Mg	Li	Zn	Cu	Mn	Sn	Ti	Zr	Ni	Pd
Basicity of M-R		Strong	Strong	Weak	Weak	Very weak		Weak	Very weak		
Substrate	Product										
R'COCl	R ₂ R'COH	(+)	(+)								
	RCOR'			++	+	+	+				
R'CHO	RR'CHOH	+	+	+	(+)	-		+	+		
Epoxide	R(CH ₂) ₂ OH	+	+		+	-					
R'I	R-R'	[+]	(+)		+	(+)	-	+			
R'COR''	RR'R''COH	+	+	+	-	-		(+)	(+)		
R'COOR''	R ₂ R'COH	+	+	-	-						
CO ₂	RCOOH	+	+	-	-						
R'COO ⁻	RCOR'	+	+	-	-						
R'C=NR''	RR'CNHR''	+	+								
R'CN	RCOR'	+	(+)								

$R'CONR''_2$	RCOR'	+		-	-						
$R'C(OR'')_3$	RCOR'	+		-	-						
C=CHal	RC=C	-	-	-	+						
Regiosel. by conjugated add.		C-1 (C-3)	C-1	C-1	C-3						
Stereosel.								Erythro	Erythro		

- ++ Best choice
- + Reacts. The reaction is suitable for preparative purposes
- (+) Reacts sluggishly or give modest yields
- [+] Reaction only suitable in certain cases
- Does not react

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Description of products and Authentication of New Compounds

Physical characteristics

The yield and the colour of a compound should always be given.

For crystalline compounds a melting point should always be given.

A crystalline compound should be recrystallized until the m.p. is constant. The recrystallization medium should always be given.

Additional physical characteristics are NMR IR UV and MS-spectra. Chiral compounds can be described by their optical rotation.

New Compounds

It is the responsibility of authors to provide fully convincing evidence for the homogeneity and identity of all compounds they claim as new.

For new compounds a correct elemental analysis is required. An accurate mass measurement of a molecular ion does not provide evidence of the purity of a compound and must be accompanied by independent evidence of the homogeneity.

Where elemental analytical data are not available, appropriate evidence which is convincing to an expert in the field may be acceptable, but authors should include, for the referees, a brief explanation of the special nature of their problem.

A compound is considered as new if

- It has not been prepared before
- If it has been prepared before but not adequately purified
- If it has been purified but not adequately characterized
- If, earlier, it has been assigned an erroneous constitution

If a compound has been described before the reported yield and m.p. and the reference should be given.

Physical Data for Journals

Characterization of a new compound in order of importance (order and format used in journals shown below):

m. p. or b. p.
 R_f (solvent)
¹H-NMR (xx MHz, solvent)
 Elemental analysis
¹³C-NMR
 MS
 [α]_D
 IR
 UV
 MS

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Litteratur

General Organic Chemistry

J. March, Advanced Organic Chemistry, Wiley, New York 1985.

nn, Functional Group Interconversion

F. A. Carey and R. I. Sundberg, Organische Chemie, Verlag Chemie, Weinheim 1995.

Smith xx

R. C. Larock, Comprehensive Organic Transformations, Verlag Chemie, Weinheim 1989.

S. R. Sandler and W. Karo, Organic Functional Group Preparations Vol. 1-3, Academic Press, New York 1986.

I. T. Harrison and S. Harrison, Compendium of Organic Synthetic Methods, Vol. 1-3, Wiley, New York 1977.
nn, Umwandlung Organischer Verbindungen.

Oxidation-Reduction

E. H. Smith, Comparative Oxidation of Organic Compounds, Wiley, New York 1989.

M. Hudlicky, Reductions in Organic Chemistry, Wiley, New York 1986.

Metallorganic Chemistry

Robert H. Crabtree, The Organic Chemistry of the Transition Metals, Wiley, New York 1988.

L. Hegedeus, xx

C. Elschenbroich and A. Salzer, Organometallics, Verlag Chemie Weinheim 1989

M. Schlosser, Organometallics in Synthesis, Wiley, New York 1994

Protective Group Chemistry

T. W. Greene, Protective Groups in Organic Synthesis, Wiley, New York

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Solvents

- C. Reichardt, *Solvents and Solvent Effects in Organic Chemistry*, 2'd Edition, Verlag Chemie, Weinheim 1988.

Technics in Synthesis

- D. D. Perrin and W. L. F. Amarego, *Purification of Laboratory Chemicals*, 3' Edltion, Pergamon, Oxford 1987.
L. Brandsma and H. Verkruijsse, *Preparative Polar Organometallic Chemistry 1*, Springer, Berlin 1987.
L. Brandsma and H. Verkruijsse, *Preparative Polar Organometallic Chemistry 2*, Springer, Berlin 1990.
C. F. Lane and G. W. Kramer, *Handling Air-sensitive Reagents*, *Aldrichimica Acta*, (1977) 11.
G. B. Gill and D. A. Whitting, *Aldrichimica Acta*, 19 (1986) 31
J. Leonard, B. Lygo and G. Proctor, *Advanced Practical Organic Chemistry*, Blackie Academic and Professional, London 1995.
M. Schlosser, *Organometallics in Synthesis A Manual*, Wiley, New York 1994
D. W. Mayo, R. M. Pike and S. S. Butcher, *Microscale Organic Laboratory*, Wiley, New York

Safety

- G. Lunn and E. B. Sansone, *Destruction of Hazardous Chemicals in the Laboratory*, Wiley, New York 1994.

Writing

- Howard M. Kanare, *Writing the Laboratory Notebook*, American Chemical Society, Washington 1985.
Robert Schoenfeld, *The Chemist's English*, Verlag Chemie, Weinheim 1990.
J. S. Dodd, ed., *The ACS Style Guide, A Manual for Authors and Editors*, American Chemical Society, Washington 1986

58 (E)-A,B-UNSATURATED ALDEHYDES**P**

Starting mat. Zinc, 1,2-Dibromoethane, t-Butylamine, Pentane, n-Butyllithium, Diisopropylamine, Tetrahydrofuran, Celit, Chlorotrimethylsilane, Zinc chloride (anhydrous), Diethyl ether

Equipment

Scale Aftales

Description *Stereoselective synthesis of E-a,b unsaturated aldehydes using Peterson olefination***92 2,2-DIMETHYL-1-PHENYL-1-PROPANONE****X**

Starting mat. n-Butyllithium, Thiophenol, Copper(I) iodide, t-Butyllithium (1.7 M in Tetrahydrofuran), Benzoyl chloride, Sodium thiosulfate, Hydrochloric acid (0.1 M)

Equipment Soxhlet (til CuI)

Scale 25 %

Description Copper catalyzed 1,4-addition

98 CIS-2-ETHYL-1-HYDROXY-2-PENTENE**S**

Starting mat. Diisobutylaluminium hydride, 1,2-Dimethoxyethane (Monoglym), Methylithium (1.4 M in Diethyl ether), 3-Hexyne, Paraformaldehyde

Equipment

Scale 10 %

Description The synthesis constitutes a general and highly flexible method for the preparation of substituted alkenes from alkynes by hydroalumination, metalation followed by electrophilic addition.

The synthesis exemplifies the use of DIBAL and methylithium under inert conditions.

99 TRANS-2-HEPTENE ACID**S**

Starting mat. Diisobutylaluminium hydride, Methylithium (1.4 M in Diethyl ether), Heptane, 1-Hexyne, Carbon dioxide

Equipment

Scale 20 %

Description The synthesis constitutes a general and highly flexible method for the preparation of substituted alkenes from alkynes by hydroalumination, metalation followed by electrophilic addition.

The synthesis exemplifies the use of DIBAL and methylithium under inert conditions.

138 ETHYL ACETOACETAT ETHYLENE ACETAL**M**

Starting mat. Ethyl acetoacetate, Ethylene glycol, 4-Toluene sulfonic acid, Toluene

Equipment Mikrosyntesesæt

Scale 100 %

Description Microsynthesis

139 TRIPHENYLMETHANOL**M**

Starting mat. Bromobenzene, Magnesium, Iodine, Benzophenone

Equipment Mikrosyntesesæt, Craigrør

Scale 100 %

Description Microsynthesis

141 2-N-BUTYL-2-METHYLCYCLOHEXANONE**S**

Starting mat. Sodium, 2-Methylcyclohexanone, Ethyl formate, Potassium, Ammonia (gas), Iron(III) nitrate nonahydrate, 1-Bromobutane

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Scale 5 %

Description The synthesis constitutes a general method for the regioselective alkylation of ketones in the most substituted alpha-position, not using thermodynamic control but using the introduction of an assisting formyl group. Subsequent di-deprotonation produces a stable and an unstable carbanionic centre, the latter being the most reactive. Alkylation and de-formylation then affords the desired pure regioisomer.

The synthesis exemplifies handling of sodium and the even more reactive potassium, the handling of liquid a

142 ETHYL 6-METHYLPYRIDIN-2-ACETATE**S**

Starting mat. 2,6-Dimethylpyridine, Potassium, Ammonia (gas), Diethyl carbonate
Equipment Omrører til ammoniakreaktionen, CO₂-Svaler
Scale 5 %
Description The synthesis illustrates handling of low boiling liquids and potassium

150 (R,R)-1,2-DIAMINODIPHENYLETHANE**S**

Starting mat. Benzil, Acetic acid, Ammonium acetate, Cyclohexanone, Argon, Tetrahydrofuran, Ammonia (gas), Lithium, Ammonium chloride, Tartaric acid (L)-(+)
Equipment Mekanisk omrører, CO₂-svaler
Scale 9 g Spiroforbindelse
Description The compound prepared can be used for the preparation of 1,2-diaminodiphenylethane which in its turn is used for the preparation of chiral ligands of the salene family for asymmetric synthesis. The chiral ligands were introduced by Jacobson and Kazuki. Jacobsen presents a Lind-Larsen lecture in Copenhagen in January 2004. The current method for the resolution of the two enantiomers of 1,2-Diaminodiphenylethane has been published by a Danish inorganic chemist from the University of Copenhagen.

184 MITSONOBU INVERSION OF STERICALLY HINDERED ALCOHOLS**P**

Starting mat. Menthol, 4-Nitrobenzoic acid, Triphenylphosphine, Diethyl azodicarboxylate, Tetrahydrofuran
Equipment Chromatografisøjle
Scale 100 %
Description Inversion of OH by Mitsunobu reaction

186 ISOPINOCAMPHEOL 2**P**

Starting mat. α -Pinene, Tetrahydrofuran, Borane (1M in Tetrahydrofuran), Sodium perborate, Nitrogen
Equipment Kuglerørsdestillationsapparat
Scale 20 %
Description Preparation and application of a versatile chiral boron hydride

243 1-BENZOYL-1,2-DIMETHYL-BUTAN-2-OL, PREPARATION OF QUATERNARY CE**T**

Starting mat. α -Bromopropiophenone, 2-Butanone, Tetrahydrofuran, Chromium(II) chloride, Lithium iodide, Argon, Sodium chloride, Ammonium chloride
Equipment
Scale 20 mmol
Description The chromium based Reformatzky reaction

252 A-AMINO ACIDS**S**

Starting mat. Methyl glyoxylate, Benzyloxyamine (sd), Nitrogen, Alkyl iodide, Triethylborane
Equipment
Scale Aftales
Description Radical reaction with triethylborane

268 (S)-1-PHENYL-1-PROPANOL (STEREOSELECTIVE ADDITION TO ALDEHYDES)**S**

Starting mat. TADDOL (sd), Titanium(IV) chloride, Tetrahydrofuran, Titanium(IV) isopropoxide, Benzaldehyde (or another aldehyde), Diethylzinc (or another dialkylzinc), Ammonium chloride, Celit
Equipment Chromatografisøjle, Polarimeter
Scale 100 %
Description Stereoselective addition of dialkyl zinc to aldehydes mediated by titaniumalkoxide

299 3-BROMO-4-(4-ETHOXYCARBONYLBENZOYL)THIOPHENE**X**

Starting mat. Lithium, Naphthalene, Manganese(II) chloride (anhydrous), Tetrahydrofuran, 3,4-Dibromothiophene, 1,2-Dibromoethane, Ethyl 4-iodobenzoate, Tetrakis(triphenylphosphine) palladium, Sodium thiosulfate

Equipment

Scale 20 %

Description The preparation and use of Rieke manganese for chemoselective acylation of aromatics

318 1-PHENYLIMIDAZOL**X**

Starting mat. Imidazole, 4-Methylphenylboronic acid, Copper(I) chloride, Tetramethylethylenediamine

Equipment

Scale 5 mmol

Description N-Arylation

324 1-(3-BROMO-4-METHOXY)-BENZYL-PYRAZOLE**S**

Starting mat. 1-(4-Methoxybenzyl)pyrazole, Bromine, Sulfur dioxide

Equipment CO2 Svaler

Scale 1 g 1-(4-Methoxy)-benzyl-pyrazol

Description The synthesis constitutes a highly regioselective bromination using bromine dissolved in liquid SO₂.The synthesis exemplifies the handling of liquid SO₂.**329 TERT-BUTYL PHENYLACETATE 1****S**

Starting mat. n-Butyllithium, Nickel bromide, Tetrahydrofuran, Iodobenzene, t-Butyl acetate, Diisopropylamine

Equipment

Scale 100 %

Description Nickel mediated synthesis of tertiary butyl esters

339 1-(4-BROMOPHENYL)PROPANOL**S**

Starting mat. Cerium(III) chloride, Diethylzinc, 4-Bromobenzaldehyde, Chlorotrimethylsilane, Toluene

Equipment Schlenk-udstyr, Handskeboks, Doseringspumpe

Scale 1.5 mmol

Description Cerium catalyzed nucleophilic addition. Handling of diethylzinc and cerium(III) chloride

347 CYCLOPROPANES A**S**

Starting mat. 1-Hexyne, Toluene, Diisobutylaluminium hydride in toluene, Diethylzinc in toluene, Diiodomethane, Bromine, Tetrahydrofuran, Palladium(II) acetate, Triphenylphosphine, Zinc chloride (anhydrous), 4-Chloriodobenzene

Equipment Rysteapparat med 24 brønde

Scale 20 %

Description The synthesis constitutes a general and highly flexible method for the preparation of substituted cyclopropanes. Alkynes used as starting materials are first hydroaluminated. The vinylalanes thus formed are then cyclopropylated by the Simmons-Smith procedure to give cyclopropylalanes. In these the aluminium can be replaced with electrophiles or cross-coupled under transition metal catalysis. In order to get proper reactivity and selectivity the aluminium can be replaced with other metals before electrophilic addition or cross-coupling.

352 PHENYLPHENANTHRIDINE**S**

Starting mat. t-Butyllithium (1.7 M in Tetrahydrofuran), Fluorobenzene, Benzonitrile (or another nitrile), Tetrahydrofuran

Equipment Digitalt lavtemperaturtermometer

Scale 100 %

Description The phenanthridine skeleton is present in several biologically active compounds. The present synthesis constitutes a new, one-pot procedure giving access to a series of substituted phenanthridines (You can try to unravel the mechanism). The protocol is the most simple of all published methods.

The synthesis exemplifies the handling of tert butyllithium, careful temperature and time control and reaction in an inert atmosphere.

359 (3S,4R)-3-HYDROXY-4-(4-METHOXYPHENYLAMINO)-4-(4-BROMOPHENYL)BUTAN-2-OL**S**

Starting mat. Proline (S), 4-Methoxyaniline, 4-Bromobenzaldehyde, Hydroxyacetone, Dimethylsulfoxide, Ammonium chloride

Equipment Chromatografisøjle, Polarimeter

Scale 200 %

Description Stereoselective cis dihydroxylation catalyzed by osmium after Sharpless

372 SILYL ENOLATES (KINETIC AND THERMODYNAMIC)**S**

Starting mat. 2-Methylcyclohexanone, Propiophenone, Tetrahydrofuran, Diisopropylamine, n-Butyllithium, Chlorotrimethylsilane, Sodium hydrogencarbonate, Pentane

Equipment NMR

Scale 5 mmol ketone

Description The handling of kinetic and thermodynamic silylenolates

401 1-HYDROXYMETHYL-4-(1-METHYLCYCLOPROPYL)-1-CYCLOHEXENE**X**

Starting mat. (S) (-)-Perillyl alcohol, Tri-isobutylaluminium (hos MB), Argon, Diiodomethane

Equipment Kuglerørsdestillationsapparat

Scale Aftales

Description The reaction is a selective cyclopropanation of an unactivated alkene using trialkylaluminium.

The method illustrates the handling of highly pyrophoric triisobutylaluminium.

407 N-ARYLATION**P**

Starting mat. Sodium hydride, Pinacol, Argon, Phenylcarbonimidic dichloride (Phenylisocyanide dichloride), Tetrahydrofuran, Triethylamine, Hydroxylammonium chloride, Benzenesulfonyl chloride, Arylmagnesium chloride

Equipment

Scale 25 %

Description A general procedure for Pd catalyzed N-arylation

409 TRIPHENYLPHOSPHINE COPPER HYDRIDE (STRYKERS REAGENT) B**S**

Starting mat. Copper(II) acetate, Triphenylphosphine, Benzene, Diphenylsilane, Acetonitrile

Equipment Schlenk, Schlenk-filter

Scale 100 %

Description Preparation of a versatile but sensitive reducing agent which can be used for 1.3-reduction.

410 VINYLBORONIC ESTER**P**

Starting mat. Sodium borohydride, 2,5-Dimethylhexadiene, Diglyme, Dimethyl sulphate, Phenylacetylene (or 4-Tolylacetylene, Formaldehyde (37 % aq), Diethanolamine

Equipment

Scale Aftales

Preparation of a sterically hindered borane useful for preparation of vinylboranes

413 3-BROMO-2-IODOBENZONITRILE T

Starting mat. 2,2,6,6-Tetramethylpiperidine, 3-Bromobenzonitrile, n-Butyllithium, Tetrahydrofuran, t-Butyllithium (1.32 M in pentane), Zinc chloride (anhydrous), Iodine, Sodium dithionite, Ammonium chloride, Nitrogen

Equipment Kuglerørsdestillationsapparat

Scale 100 %

Description The reaction constitutes a versatile regioselective deprotonation of an aromatic compound followed by addition of an electrophile to give a multifunctionalized aromatic compound.

The procedure illustrates the handling and distillation of highly pyrophoric dialkylzinc. The method also illustrates the handling of n-butyllithium and LiTMP.

414 N-(2-DICYCLOHEXYLPHOSPHINOPHENYL)-2,2,6,6-TETRAMETHYLPYPERIDINE S

Starting mat. n-Butyllithium, 2,2,6,6-Tetramethylpiperidine, Iodobenzene, t-Butyllithium, Tetrahydrofuran, Dicyclohexylchlorophosphine, Nitrogen

Equipment

Scale Aftales

Description Preparation of an efficient ligand for palladium

422 2-METHYL-3-PHENYLPROPANOIC ACID S

Starting mat. n-Butyllithium, Diisopropylamine, Tetrahydrofuran, Nitrogen, 3-Phenylpropanoic acid, Iodomethane

Equipment

Scale 33 %

Description Dimetalation of a carboxylic acid followed by alpha-alkylation

425 METHYL 4-HEXYLBENZOATE S

Starting mat. Argon, Methyl 4-chlorobenzoate, Iron(III) acetylacetonate, Tetrahydrofuran, N-Methylpyrrolidone, Hexylmagnesiumbromide (2 M in diethyl ether)

Equipment Schlenk

Scale 100 %

Description A Corriu-Tamao type cross-coupling.

428 1-BENZYLOXY-5-TRIMETHYLSILYLETHYNYLPYRAZOLE S

Starting mat. 1-Benzyloxy-5-iodopyrazole (sd), Trimethylsilylacetylene, Tetrahydrofuran, Nitrogen, Triethylamine, Palladium bis(triphenylphosphine)dichloride, Triphenylphosphine, Copper(I) iodide

Equipment

Scale Aftales

Description A Sonogashira cross-coupling.

432 CYCLOPROPANATION S

Starting mat. Cinnamyl alcohol, t-Butyldimethylsilylchloride, Triethylamine, Diethylzinc (1.0 M in hexanes), Nitrogen, Trifluoroacetic acid, Calcium hydride, Diiodomethane

Equipment Søjle

Scale Aftales

Description A versatile cyclopropanation method.

442 [1-(S)-BENZYL-5-(4-(R)-BENZYL-2-OXO-OXAZOLIDIN-3-YL)-4-(S)-HYDROXYMET S

Starting mat. Diisopropylamine, n-Butyllithium, Tetrahydrofuran, 3-Acetyl-4-(R)-benzyl-oxazolidin-2-one (see 157 A), Cbz, PMB-Phe-OH (see 439), N,N-carbonyldiimidazole, Diethyl zinc in hexane, Diiodomethane, Paraformaldehyde, Ammonium chloride

Equipment Doseringspumpe

Scale Aftales

A ketomethylene peptide mimetic.

454 DIMETHYLITANOCENE, TEBBE REACTION**S**

Starting mat. Nitrogen, Titanocene dichloride, Toluene, Methylmagnesium chloride, Tetrahydrofuran, Ammonium chloride

Equipment

Scale 10 %

Description Preparation of dimethyltitanocene, an excellent substitute for the Tebbe reagent

470 1,2,3,4-TETRAFLUOROBENZEN-5,6-DISULFONIC ACID MONOHYDRATE**S**

Starting mat. Nitrogen, 1,2-dibromotetrafluorobenzene, n-Butyllithium, Sulphur dioxide, Hydrogenperoxide aq 25 %, Amberlite IR 120, Potassium hydroxide

Equipment Schlenk

Scale Aftales

Description Generation of an aromatic dianion followed by reaction with a (gaseous) electrophile

486 FMOC-ARG(BOC)2-OH**P**

Starting mat. Brian Lohse

Equipment

Scale Aftales

Description Otherwise difficult protection of arginine for subsequent peptide synthesis is effected using ornithine as the starting material

487 FMOC-ARG(BOC)2-OH B**P**

Starting mat. Brian Lohse

Equipment

Scale Aftales

Description Otherwise difficult protection of arginine for subsequent peptide synthesis is effected using ornithine as the starting material

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2009 ©